


Belsorp Adsorption/Desorption Data Analysis Software - Ver 6.1.0.4

BEL Japan, Inc.



BJH-Plot

Adsorption branch

Adsorptive N2

Adsorption temperature 77 [K]

○ 1207682-GMSI-N-2-N2 Test.DAT

1207682-GMSI-N-2-N2 Test
Catalyst Evaluation Department
Degassed @305C for 3h
Vacuum degree before measurement: 7.161E-5Pa

Sample weight 0.3142 [g]

Date of measurement 12/02/27

Saturated vapor pressure 88.493 [kPa]

Time of measurement 25:24:49

V_p 0.4292 [cm³ g⁻¹]

a_p 190.89 [m² g⁻¹]

احمد فضل الله